



29 CFR 1910.1030 – Bloodborne Pathogens

Universal Precautions is an approach to infection control. According to the concept of Universal Precautions, all human blood and certain human body fluids are treated as if known to be infectious for HIV, HBV, and other bloodborne pathogens.

Five Steps to Efficiently Contain and Remove a Bloodborne Pathogens Spill with the Safetec® Universal Precaution Kit*

*All personnel must be trained in the use of the spill kit before any cleanup process is performed

(To help remember the steps, think of the acronym **PCRDDS**)

1. **PREVENT** - to prevent direct contact with the spill during the cleanup process and minimize the risk for infection, disposable protective clothing must be worn. This includes *gloves, combo mask/safety shield and protective apron*.
2. **CONTAIN & REMOVE** - to efficiently contain and remove the spill, a *solidifying agent (Safetec® Red Z®)* and *scoop/scrapper* are to be used and the spill contents are to be disposed in a *red biohazard bag*.
3. **DISINFECT** - to properly disinfect the contaminated surface area, use the (SaniZide Plus®) germicidal wipe. (This process is critical for preventing any further contamination)
4. **DISPOSE** - discard the germicidal wipe, gloves, safety shield and apron into *red biohazard bag* and dispose of bag as required by your local, state and federal regulations.
5. **SANITIZE** - use the *p.a.w.s® wipe* to sanitize hands and use soap and water after immediate cleanup is done

For a more detailed description of this process, please view the **Safetec Video Demonstration** at <http://safetec.com/distributor-resources/> or scan the QR code to the right.



Scan the QR code to learn how to use the Universal Precaution Kit

Helpful Tips to Remain OSHA Compliant

- Think safety culture – not just compliance
- Compliance is the key to prevention
- Do not guess the information or make assumptions
- Plan for an inspection ahead of time
(Do not be in a position of figuring it out as you go)

*For more information about OSHA and helpful tips, visit www.safetec.com/OSHA-Compliance



References

Occupational Safety and Health Administration. Bloodborne Pathogens and Needlestick Prevention. Accessed March 2014 from: <http://www.osha.gov/SLTC/bloodbornepathogens/index.html>

Occupational Safety and Health Administration. (2011) OSHA Fact Sheet: OSHA's Bloodborne Pathogens Standard. Accessed March 2014 from: https://www.osha.gov/OshDoc/data_BloodborneFacts/bbfact01.pdf

Ippolito, D. & Sierra G. (2014, March). How to Prepare for & Manage an Aggressive OSHA Inspection [Webinar].

