



Fast-Acting



Self-Dissolving



No Splashing

## Drop-In Pacs

### Advanced Self-Dissolving Paper Technology

Drop-In Pacs are a solidifier with an advanced self-dissolving technology to instantly absorb and solidify fluid waste. Containing the same effective super absorbent polymers as used in Red Z® and Green-Z® bottles, these slender pouches are made with a fast-acting innovative water-soluble paper.

#### Ideal for use in:

- Kick Buckets
- Operating Rooms
- Biohazard Bags
- Sponge Receptacles
- Substance Transport
- Suction Canisters
- Wound Care Therapy

### Directions for Use

**Suction Canisters:** Insert entire Drop-In Pac through port opening of suction canister (before or after surgical procedure). Allow sufficient time for encapsulation before handling. Dispose of infectious waste\* or chemical spill\*\*\*.

**Biohazard Bags or Kick Buckets:** Place drop-in pac in the bottom of a biohazard bag or other container to control fluid before or after spills occur. Dispose of infectious waste\* or chemical spill\*\*\*. Clean and disinfect contaminated areas after use.



**OSHA Compliance 1910.1030 – Bloodborne Pathogens Standard**

Universal Precautions is an approach to infection control. According to the concept of Universal Precautions, all human blood and certain human body fluids are treated as if known to be infectious.

Order No.	Description	Contains Chlorine	Case Qty.	Time to Open	Time to Absorb	Capacity*
41108	Red Z® Drop-In (up to 750cc)	✓	300/case	10-25 sec	1-3 min	750-900cc
41175	Red Z® Drop-In (up to 1500cc)	✓	100/case	10-25 sec	2-6 min	1500-1800cc
42008	Green-Z® Drop-In (up to 750cc)	-	300/case	10-25 sec	1-3 min	750-900cc
42175	Green-Z® Drop-In (up to 1500cc)	-	100/case	10-25 sec	2-6 min	1500-1800cc

# SOLIDIFY FLUID WASTE QUICKLY & SAFELY

*Self-Dissolving paper technology is ideal for improving fast-paced operating room turnover and worker safety*

*Pouch is specifically designed to fit in the porthole of suction canisters*



\*Data is an approximation of values. Blood and urine solidification are based on the industry standard 0.9% NaCl saline solution.

\* Dispose of infectious and hazardous waste in accordance with all local, state and federal regulations.

\*\*\* Careful consideration should be given to chemical spill cleanup, i.e., room size, temperature, quantity and concentration(s) of the solution(s).



**Made in the USA**

**Safetec of America, Inc.**  
1-800-456-7077 | [safetec.com](http://safetec.com)