



Fast-Acting



Chlorine-Free



No Splashing

## Green-Z®

### Fluid Solidifier



At times, workers may be unintentionally exposed to potentially infectious blood and body fluid spills along with other harmful substances such as toxic chemicals, mercury spills, and chemotherapy agents. To help reduce those exposure risks, rely on Green-Z® for improved safety with spill cleanup. Green-Z® is a chlorine-free powder fluid solidifier that transforms liquids into a gel-like mass, providing a reliable and efficient solution that is safe for incineration. Green-Z® reduces the risks associated with improper spill cleanup, is user-friendly, eco-conscious, and trusted across industries.

Available in single & multi-use bottles, single-use pouches, drop-in pacs, and buckets.

- Contains super absorbent polymers (SAPs)
- Easy-to-use
- Eliminates splashing and aerosolization
- Fast-acting
- Safe for incineration
- Shelf-stable in a cool, dry place

### Ideal for use in:

- Correctional Facilities
- Healthcare & Laboratories
- Industrial & Manufacturing
- Occupational Health & Safety
- Restaurants
- Schools
- Transportation
- Anywhere spills can happen

\*Data is an approximation of values. Blood and urine solidification are based on the industry standard 0.9% NaCl saline solution. On average, 1 gram of solidifier absorbs 37.5 cc of liquid. (see charts on the following page for specific product absorbency)

# SPILL CONTROL

**Chlorine-Free For Safe Use On Chemical Spills**

**Solidify Blood & Body Fluid Spills With Ease**



\*\*\* Careful consideration should be given to chemical spill cleanup, i.e., room size, temperature, quantity and concentration(s) of the solution(s).

# Green-Z® Configurations

## Single & Multi-Use Bottles

Order No.	Size	Absorbency*	Cap	Case Qty.	Case Wt.
42003	16 oz.	Up to 22,000 cc	Needle Nose	12/case	15 lbs.
42005	2.2 oz.	3,000 cc	Standard	75/case	16 lbs.
42010	5 oz.	5,000 cc	Shaker	24/case	11 lbs.
42012	15 oz.	Up to 18,000 cc	Shaker	12/case	14 lbs.



## Single-Use & Pour-In Pouches

Order No.	Size	Absorbency*	Bags	Case Qty.	Case Wt.
42019	21 g	787 cc	-	100/case	6 lbs.



## Drop-In Pacs & Mini Pacs

Order No.	Size	Absorbency*	Time to Open	Time to Solidify	Bags	Case Qty.	Case Wt.
42004	10 g	375 cc	10-25 sec	1-3 min	100 ct.	10/case	26 lbs.
42008	20 g	750 cc	10-25 sec	1-3 min	-	300/case	16 lbs.
42029	4 g	150 cc	10-25 sec	1-3 min	500 ct.	10/case	55 lbs.
42030	0.5 g	18.75 cc	10-25 sec	1-3 min	500 ct.	3/case	5 lbs.
42031	2 g	75 cc	10-25 sec	1-3 min	-	1,000/case	7 lbs.
42175	40 g	1,500 cc	10-25 sec	2-6 min	-	100/case	12 lbs.



## Bucket

Order No.	Size	Absorbency*	Product Type	Case Qty.	Case Wt.
42220	50 lbs.	800,000 cc	Drum	1/case	55 lbs.

42220 is made to order.



Page 2  
Rev 25.10

# Green-Z® Configurations (Cont.)

## Kits Containing Green-Z®

Order No.	Description	Container	Case Qty.	Case Wt.
49725	Hazardous Drugs Spill Kit	Plastic Bag	12/case	11 lbs.
15201	Multi-Purpose Spill Kit	Poly Bag	1/case	11 lbs.
15205	Spill Leader Kit	Poly Bag	1/case	5 lbs.



## Directions for Use



### 1. Wear Personal Protective Equipment (PPE)

Wear appropriate personal protective equipment when handling any kind of spill.



### 2. Sprinkle, Pour, or Drop in

Use the Green-Z® format best suited for your spill/fluid. Allow sufficient time for solidification to occur.



### 3. Remove Solidified Spill

Dispose of your solidified spill safely. For uncontained spills, use a scoop and scraper (or similar tool) and place in an appropriate container.



### 4. Clean and Disinfect

Clean and disinfect contaminated areas after use and wash thoroughly.



### OSHA Compliance 1910.1030 – Bloodborne Pathogens Standard

Universal Precautions is an approach to infection control. According to the concept of Universal Precautions, all human blood and certain human body fluids are treated as if known to be infectious.

### Expiration Dating

Green-Z® is shelf stable and does not require an expiration date. In order for the product to remain shelf stable, it must be kept in a cool, dry place and in its original packaging. If contents are hard or compact, simply shake the product to loosen, before use.

*\*Dispose of infectious and hazardous waste in accordance with all local, state and federal regulations.*