

Cleaner

Deodorizer

Disinfectant

SaniZide Plus® Germicidal Solution

Kills 99% of Bacteria #



SaniZide Plus® is a one-step*¹ germicidal cleaner and deodorant designed for general cleaning, disinfecting and controlling non-pathogenic mold and mildew on hard, non-porous inanimate surfaces. Quickly removes dirt, grime, non-pathogenic mold and mildew, food residue, blood and other organic matter commonly found in health care facilities. It also eliminates odors leaving hard non-porous surfaces smelling clean and fresh. Use where odors are a problem. This product is effective against 48 microorganisms on hard non-porous surfaces. (See table for exact kill times.)

Intended for use in:

- Ambulances
- Ambulatory Surgical Centers
- Animal Care Facilities
- Bathrooms
- Hospitals
- Laboratories
- Medical or Dental Offices
- Nursing Homes
- Office Buildings
- Other Healthcare Related Institutions
- Salons
- Veterinary Clinics

Order No.	Size	Product Type	Carton	Case Qty.	Case Wt.
34800	4 oz.	Spray Bottle	-	24/case	8 lbs.
34805	16 oz.	Spray Bottle	-	12/case	15 lbs.
34810	32 oz.	Spray Bottle	-	6/case	15 lbs.
34815	1 Gallon	Refill Bottle	-	4/case	38 lbs.
34855	55 Gallon	Drum	-	1/case	520 lbs.

CLEAN AND DEODORIZE IN ONE EASY STEP*¹

**Effective against COVID-19 and Hantavirus
on hard non-porous surfaces**

**Meets OSHA Bloodborne Pathogen
Standard for HIV, HBV, and HCV**

¹ "One Step" disinfectant claim is limited to association with non-food contact surface use.

* when used according to disinfection directions



SaniZide Plus® Germicidal Solution Efficacy Sheet

May be used on hard, non-porous surfaces in:

Hospitals and healthcare settings, including: operating rooms, nursing homes, ICU areas, and patient rooms. **Animal Housing/Holding Areas:** crates, cages, kennels, animal life science laboratories, grooming facilities, pet stores, and veterinary clinics. **Other Use sites:** athlete/recreational facilities, homes, lodging establishments, schools, food processing plants, shower rooms, bathrooms, retail businesses, laminated surfaces, empty garbage bins, tubs, urinals, and walls.

Virucidal Activity (Viruses):

Test Organism	Contact Time
Virus Enveloped:	
HIV-1 (associated with AIDS).....	1 minute
Avian Flu A (Strains H9N2 and H3N2).....	2 minutes
Hepatitis B Virus (HBV).....	1 minute
Hepatitis C Virus (HCV).....	1 minute
Duck Hepatitis B Virus (DHBV).....	1 minute
Bovine Viral Diarrhea Virus (BVDV).....	1 minute
Human Coronavirus.....	2 minutes
SARS Associated Coronavirus (ZeptoMetrix).....	2 minutes
SARS-Related Coronavirus 2 [SARS-CoV-2] USA-WA1/2020 BEI NR-52281.....	1 minute
Rabies Virus.....	30 seconds
Paramyxovirus (Mumps).....	3 minutes
Virus Large Non-Enveloped:	
Rotavirus.....	3 minutes
Virus Small Non-Enveloped:	
Norovirus (Norwalk virus).....	30 seconds
Feline Calicivirus.....	30 seconds
Canine Parvovirus.....	10 minutes
Poliovirus Type 1.....	10 minutes
Rhinovirus type 39.....	3 minutes
Hepatitis A Virus (HAV) (University of Ottawa).....	10 minutes

Bactericidal Activity (Bacteria):

Test Organism	Contact Time
Pseudomonas aeruginosa.....	3 minutes
Staphylococcus aureus.....	3 minutes
Salmonella (choleraesuis) enterica (ATC 10708).....	3 minutes
Escherichia coli (E. coli).....	3 minutes
Escherichia coli O157:H7.....	3 minutes
Streptococcus pyogenes (Necrotizing Fasciitis-Group A) (V.A. Medical Center Isolate 04001).....	3 minutes
Listeria monocytogenes.....	3 minutes
Yersinia enterocolitica.....	3 minutes
Enterococcus faecium.....	3 minutes
Corynebacterium ammoniagenes.....	3 minutes
Salmonella (typhi) enterica (ATCC 6539).....	3 minutes

Antibiotic-Resistant (Strains of) (Antibiotic-Resistant Bactericidal Activity):

Test Organism	Contact Time
Methicillin resistant Staphylococcus aureus (MRSA).....	3 minutes
Vancomycin resistant Enterococcus faecalis (VRE).....	3 minutes
Vancomycin intermediate resistant Staphylococcus aureus (VISA).....	3 minutes
Methicillin resistant Staphylococcus epidermidis (MRSE).....	3 minutes
Community Associated Methicillin resistant Staphylococcus aureus (CA-MRSA) (NRS 123) Genotype USA400.....	3 minutes
Community Associated Methicillin resistant Staphylococcus aureus (CA-MRSA) (NRS 384) Genotype USA300.....	3 minutes

Tuberculocidal

Test Organism	Contact Time
Mycobacterium tuberculosis BCG (Mycobacterium bovis).....	5 minutes

Mildew/ Fungi

Test Organism	Contact Time
Trichophyton interdigitale, formerly Trichophyton mentagrophytes (athlete's foot fungus).....	10 minutes

Safetec's SaniZide Plus Spray/Liquid (EPA Reg # 1839-83-67161) has demonstrated effectiveness against viruses similar to Hantavirus on hard, non-porous surfaces. Therefore, this product can be used against Hantavirus when used in accordance with the directions for use against Human Rotavirus or Rhinovirus Type 39 on hard, non-porous surfaces. Refer to the CDC website at <https://www.cdc.gov/hantavirus/hcp/clinical-overview/hps.html> for additional information.