



Animal Care

- Removes dirt, grime, blood, urine, and fecal matter
- Reduces odor-causing bacteria
- 99.9% effective in killing microorganisms and germs



Safetec SaniZide Pro 1[®] 1-Minute Surface Disinfectant

Protect animals from harmful germs and make patient turnover times even faster with Safetec’s SaniZide Pro 1[®] Surface Disinfectant. This ready-to-use, alcohol-based formula combines cleaning and disinfecting in one simple step with true, one-minute kill claims that are effective against 48 microorganisms, including **Canine Parvovirus, Canine Distemper virus, Enterococcus, Monkeypox, Streptococcus, & Trichophyton (Ringworm)**.

SaniZide Pro 1[®] is non-abrasive and safe for repeated use on hard, non-porous surfaces associated with non-food contact countertops, exam tables, cages, kennels, animal equipment, devices found in animal housing facilities, and more.

Ideal for use in:

- Animal Hospitals
- Animal Life Science Laboratories
- Animal Transportation Vehicles
- Breeding Establishments
- Grooming Establishments
- Kennels
- Pet Shops
- Tack Shops
- Veterinary Clinics
- Zoos


Order No.	Size	Product Type	Carton/Canister	Case Qty.	Case Wt.
35910	32 oz.	Spray Bottle	-	6/case	14 lbs.
35923	6" x 6.75"	Wipe Canister	160 ct.	12/case	22 lbs.
35926	8" x 8"	Individual Wipe	50 ct.	6/case	8 lbs.
35927	8" x 8"	Individual Wipe	Bulk	500/case	15 lbs.
35941	5 Gallon	Carboy	-	1/case	40 lbs.
35955	55 Gallon	Drum	-	1/case	520 lbs.

ANSI Supplemental Product

EPA Reg. No.'s: 88494-3-67161 (spray)
88494-4-67161 (wipes)

SaniZide Pro 1[®] contains Ethanol (CAS 64-17-5) as a disinfectant, which in California is considered to be carcinogenic to humans when ingested in the form of an alcoholic beverage [IARC 1]. SaniZide Pro 1[®] also contains the following ingredients which do not appear on any listing of chemicals of concern as designated by the State of California or the European Union: Purified Water (CAS 7732-18-5) as a diluent, Propylene Glycol (CAS 57-55-6) as a solvent, Citric Acid (CAS 77-92-9) as a pH adjuster, Decyl Glucoside (CAS 68515-73-1) as a surfactant, Didecylidmonium Chloride (CAS 7173-51-5) as a disinfectant, Lauryl Glucoside (CAS 110615-47-9) as a surfactant.

Page 1
Rev 26.03



Animal Care Continued

Directions For Use

Veterinary Clinics, Animal Life Science Laboratories, Animal Care/Housing: Remove all animals and feed from premises, vehicles, and enclosures. Remove all litter and manure from floors, walls, and surfaces of facilities and fixtures occupied or traversed by animals. Empty all troughs, racks and other feeding and watering appliances.

Thoroughly clean all surfaces with soap or detergent and rinse with water. Using either SaniZide Pro 1[®] Spray or Wipes, wet surfaces for 1 minute. Immerse all handling and restraining equipment as well as forks, shovels, and scrapers used for removing litter and manure.

Ventilate buildings, vehicles, and other closed spaces. Do not house animals or employ equipment until treatment has been absorbed, set, or dried. Thoroughly scrub treated feed troughs, racks and other feeding and watering appliances with soap or detergent, then rinse with potable water before reuse.

Nonrefillable container. Do not reuse or refill this container. Offer for recycling if available, or discard empty container in trash.



SaniZide Pro 1[®] Efficacy Information

May be used on hard, non-porous surfaces in Animal Housing/Holding Areas: Animal Hospitals, Animal Life Science Laboratories, Animal Transportation Vehicles, Breeding Establishments, Grooming Establishments, Kennels, Pet Shops, Tack Shops, Veterinary Clinics, Veterinary Examination Rooms, Veterinary Operating Rooms, Veterinary Waiting Rooms, Veterinary X-Ray Rooms, Zoos.

Test Organism

†Avian Influenza A Virus	1-minute
Bordetella bronchiseptica	1-minute
Bovine Diarrhea Virus	1 minute
† Canine distemper virus	1-minute
†† Canine parvovirus	1-minute
Duck Hepatitis B Virus	1-minute
Enterobacter aerogenes or Enterobacter cloacae	1-minute
Enterococcus.....	1 minute
††Enterovirus EV-D68.....	1-minute
†Feline leukemia virus	1-minute
† Monkeypox	1-minute
†SARS-CoV-2 (Coronavirus)	1-minute
Streptococcus pneumoniae	1-minute
Streptococcus pyogenes	1-minute
Trichophyton mentagrophytes (Ringworm, Athlete's Foot, Jock Itch)	1-minute
†Enveloped/lipophilic type virus ††Non-enveloped/hydrophilic type virus	

Contact Time

Safetec's SaniZide Pro 1[®] Products have demonstrated effectiveness against emerging viral pathogens similar to the Monkeypox virus on hard, nonporous surfaces when used in accordance with the directions according to the label. A list of approved products can be found on the EPA's Disinfectants for Emerging Viral Pathogens (EVPs): List Q.

Safetec's SaniZide Pro 1[®] Products have demonstrated effectiveness against SARS-CoV-2 virus on hard, nonporous surfaces when used in accordance with the directions according to the label. A list of approved products can be found on the EPA's Disinfectants for Emerging Viral Pathogens (EVPs): List N.