

DISINFECT SURFACES

Gently and Effectively

SaniZide Plus[®] Surface Disinfectant Spray

SaniZide Plus[®] is a convenient, fast-acting, multi-purpose, broad-spectrum cleaner, disinfectant, and deodorizer for environmental surfaces. Our alcohol-free formulation is a non-flammable, non-corrosive, EPA registered, quaternary ammonium compound that meets non-hazardous shipping requirements. SaniZide Plus[®] controls mold and mildew, helps prevent cross-contamination, and helps you comply with the OSHA Bloodborne Pathogens Standard, which requires the use of an intermediate level (tuberculocidal) disinfectant where gross contamination is known to have occurred.



**Effective against MRSA,
 VRE, Tuberculosis, Influenza,
 Poliovirus, and Norovirus**

Item#	Description	Case Qty.
34800	SaniZide Plus [®] 4oz. spray bottle	24/case
34805	SaniZide Plus [®] 16oz. spray bottle	12/case
34810	SaniZide Plus [®] 32oz. spray bottle	6/case
34815	SaniZide Plus [®] 1 gal. refill bottle	4/case

Ideal for use in:

- Hospitals
- Medical or Dental Offices
- Ambulances
- Ambulatory Surgical Centers
- Nursing Homes
- Laboratories
- Pharmacies
- Other Healthcare Related Institutions

SaniZide Plus® Efficacy Sheet

Virucidal Activity

Test Methods:

* U.S. E.P.A. Pesticide Assessment Guidelines, Subdivision G: Product Performance, Section 91-2(f), and Section 91-30 (d), (e), November, 1982.

† Protocols for Testing the Efficacy of Disinfectants against Hepatitis B Virus (HBV) (EPA, Federal Register, Vol. 65, No. 166, 8/25/2000, p. 51828).

‡ Protocol for Testing Disinfectants against Hepatitis C Virus using Bovine Viral Diarrhea Virus as approved by the U.S. EPA on August 15, 2002.

**Modified U.S. E.P.A. Pesticide Assessment Guidelines, Subdivision G: Product Performance, Section 91-2(f), and Section 91-30 (d), (e), November, 1982.

***Per E.P.A. Guidance Document dated 10/31/09.

Test Organism

Contact Time

*Hepatitis A Virus (HAV)	10 minutes
†Hepatitis B Virus (HBV)	5 minutes
‡Hepatitis C Virus (HCV)	5 minutes
*HIV-1 (associated with AIDS)	1 minute
*Human Coronavirus.....	2 minutes
**Norovirus (Norwalk Virus).....	30 seconds
*Poliovirus Type 1, strain Brunhilde.....	10 minutes
*SARS Associated Coronavirus (ZeptoMetrix)	2 minutes
*Avian Flu A (Strains H9N2 and H3N2)	2 minutes
***2009 H1N1 Influenza A Virus.....	2 minutes
Mycobacterium bovis BCG (Organon Teknika)	5 minutes
Trichophyton mentagrophytes	10 minutes
Staphylococcus aureus.....	3 minutes
Salmonella (choleraesuis) enterica (ATC 10708)	3 minutes
Pseudomonas aeruginosa	3 minutes
Escherichia coli (E. coli).....	3 minutes
Salmonella (typhi) enterica (ATCC 6539).....	3 minutes
Streptococcus pyogenes.....	3 minutes
Community Associated Methicillin resistant Staphylococcus aureus (CA-MRSA).....	3 minutes
Methicillin resistant Staphylococcus aureus (MRSA).....	3 minutes
Methicillin resistant Staphylococcus epidermidis (MRSE)	3 minutes
Vancomycin resistant Enterococcus faecalis (VRE).....	3 minutes
Vancomycin intermediate resistant Staphylococcus aureus (VISA).....	3 minutes

SaniZide Plus® (EPA Reg # 1839-83-67161) has demonstrated effectiveness against viruses similar to 2019 novel coronavirus (SARS-CoV-2) on hard non-porous surfaces. Therefore, this product can be used against 2019 novel coronavirus (SARS-CoV-2) when used in accordance with the directions for use against Hepatitis A virus on hard, non-porous surfaces. Refer to the CDC website (<https://www.cdc.gov/coronavirus/2019-ncov/index.html>) for additional information.

For use on:

Gurneys, isolettes, hospital beds, shower stalls, sinks, toilets, grocery carts, changing tables, cribs, trash cans, gym equipment, non-porous personal protective equipment or anywhere normally treated with a disinfectant.

Disclosure: For US State Right to Know Act, please see section 15 in our SDS.

SaniZide Plus® contains Ethanol (CAS 64-17-5) as a disinfectant, which in California is considered to be carcinogenic to humans when ingested in the form of an alcoholic beverage [IARC 1]. SaniZide Plus® also contains the following ingredients which do not appear on any listing of chemicals of concern as designated by the State of California or the European Union: Purified Water (CAS 7732-18-5) as a diluent, Propylene Glycol (CAS 57-55-6) as a solvent, Citric Acid (CAS 77-92-9) as a pH adjuster, Decyl Glucoside (CAS 68515-73-1) as a surfactant, Didecylidimonium Chloride (CAS 7173-51-5) as a disinfectant, Lauryl Glucoside (CAS 110615-47-9) as a surfactant.