



## Formaldehyde Spill Response Kit

A convenient emergency response kit for medical personnel using formalin or other formaldehyde-containing substances. It is ideal for containing small-volume spills.▪ Available in a highly visible, durable, easy-to-use hard vinyl case.

This kit includes FSC 1, a formalin spill control solidifier that aids in reducing the exposure to formaldehyde gas by immobilizing formalin and transforming it to an easily removable semi-solid mass. The resulting mass cools rapidly to minimize formaldehyde off-gassing. The capacity of FSC 1 to convert and immobilize formalin decreases as the concentration of formalin increases.

Formaldehyde Strength (% wt.)		4%	10%	15%	20%	30%	37%
<b>Absorption Capacity</b> (per 11 oz. Jar)	milliliters (mL)	1800	720	480	360	240	195
	fluid ounces (fl. oz.)	60.8	24.3	16.2	12.1	8.1	6.5

Order No.	Container	Container Dimensions	Case Qty.	Case Wt.
48625	Poly Bag	11.75" x 11.75" x 4.88"	1/case	4 lbs.
48630	Hard Case	14" x 12.625" x 5.375"	1/case	5 lbs.

### What's Included:

- 4 Clear Disposal Bags (30" x 36")\*
- FSC 1 Solidifier (11 oz. Solidifier)
- Instructions
- Open Back Gown
- Pair of Nitrile Gloves
- Respirator
- 4 Scoops/Scrapers
- 4 Twist Ties
- Vented Goggles
- 3 Wiper Pads (11" x 11")

**Please Note:** Not recommended to be used with 40% formalin or higher.

# REDUCE EXPOSURE

**Immobilize formalin and minimize off-gassing**

\* Dispose of potentially infectious materials in accordance with all local, state and federal regulations.

▪ Where a significant volume of fluid must be contained, notify emergency response personnel in accordance with facility guidelines.

**Safetec**  
of America, Inc.

**Made in the USA**

**Safetec of America, Inc.**  
1-800-456-7077 | [safetec.com](http://safetec.com)

