



# Formaldehyde Spill Response Kit

A convenient emergency response kit for medical personnel using formalin or other formaldehyde-containing substances. It is ideal for containing small-volume spills.▪ Available in a highly visible, durable, easy-to-use hard vinyl case.

This kit includes FSC 1, a formalin spill control solidifier that aids in reducing the exposure to formaldehyde gas by immobilizing formalin and transforming it to an easily removable semi-solid mass. The resulting mass cools rapidly to minimize formaldehyde off-gassing. The capacity of FSC 1 to convert and immobilize formalin decreases as the concentration of the formalin increases.

## What's Included:

- 4 Clear Disposal Bags (30" x 36")\*
- FSC 1 Solidifier (11 oz. Solidifier)
- Instructions
- Open Back Gown
- Pair of Nitrile Gloves
- Respirator
- 4 Scoops/Scrapers
- 4 Twist Ties
- Vented Goggles
- 3 Wiper Pads (11" x 11")

**Please Note:** Not recommended to be used with 40% formalin or higher.

Order No.	Container	Container Dimensions	Case Qty.	Case Wt.
48625	Poly Bag	11.75" x 11.75" x 4.88"	1/case	4 lbs.
48630	Hard Case	14" x 12.625" x 5.375"	1/case	5 lbs.

# REDUCE EXPOSURE

**Immobilize formalin and minimize off-gassing**

\*Dispose of potentially infectious materials in accordance with all local, state and federal regulations.

▪ Where a significant volume of fluid must be contained, notify emergency response personnel in accordance with facility guidelines.

