



Fast-Acting



Doesn't Expire\*\*



No Splashing

## FSC 1 Formalin Spill Control

FSC 1 is a formalin spill control solidifier that aids in reducing the exposure to formaldehyde gas by immobilizing formalin and transforming it to an easily removable semi-solid mass. The resulting mass cools rapidly to minimize formaldehyde off-gassing. It is ideal for containing small-volume spills.■

The capacity of FSC 1 to convert and immobilize formalin decreases as the concentration of formalin increases.

Formaldehyde Strength (% wt.)		4%	10%	15%	20%	30%	37%
Absorption Capacity (per 11 oz. Jar)	milliliters (mL)	1800	720	480	360	240	195
	fluid ounces (fl. oz.)	60.8	24.3	16.2	12.1	8.1	6.5

**Please Note:** Not recommended to be used with 40% formalin or higher. Careful consideration should be given to chemical spill cleanup, i.e., room size, temperature, quantity and concentration(s) of the solution(s). Dispose of infectious and hazardous waste in accordance with all local, state and federal regulations.

Order No.	Size	Product Type	Carton	Case Qty.	Case Wt.
47101	11 oz.	Shaker Top Bottle	-	12/case	11 lbs.

### Ideal for use in:

- Educational Institutions
- Funeral Homes
- Healthcare & Laboratories
- Industrial & Manufacturing
- Mortuaries
- Occupational Health & Safety
- Textile and Wood Industries

A complete solution for Formalin or Formaldehyde spills: [Safetec's Formaldehyde Spill Response Kit](#)

# REDUCE EXPOSURE

**Immobilize formalin and minimize off-gassing**



\*\*Products are shelf-stable and do not require an expiration date. For the product to remain shelf-stable, it must be kept in a cool, dry place and in its original packaging. If contents are hard/compact, shake the product to loosen before use.

■ Where a significant volume of fluid must be contained, notify emergency response personnel in accordance with facility guidelines.



Made in the USA

Safetec of America, Inc.  
1-800-456-7077 | [safetec.com](http://safetec.com)