



+



Universal Spill Response



Multi-Use

## Multi-Purpose Spill Kit

### Solidify Various Spills

This portable, easy-to-carry kit is your go-to solution for a quick and efficient response to a variety of spills. Whether dealing with chemicals, body fluids, or other hazardous substances, this kit equips you with multiple fluid solidifiers and tools for safe spill cleanup. Its compact design makes it easily accessible, ideal for a variety of environments—from laboratories to workplaces—and it conveniently fits behind a truck seat or in a spill response cart.

Meets requirements for compliance with DOT 1713 hazardous discharge; EPA 623.30 immediate action & 263.31 discharge cleanup.

#### Acid Lock Solidifier

Neutralizes and solidifies **low pH spills**: Acetic Acid, Hydrochloric Acid, Nitric Acid, and Sulfuric Acid

#### Alky Solidifier

Neutralizes and solidifies **high pH spills**: Sodium Hydroxide, Potassium Hydroxide and Ammonium Hydroxide

#### Green-Z® Solidifier

Solidifies **blood and body fluids**, and as well as **chemical\*\*\* spills**

#### Petro Lock Solidifier

Absorbs **petroleum-based spills**: Gasoline, Solvents, Crude Oil, and Paint Sludge

Order No.	Container	Container Dimensions	Case Qty.	Case Wt.
15201	Poly Bag	12" x 17.75" x 4.25"	1/case	11 lbs.

### What's Included:

- Acid Lock Solidifier (15 oz. Bottle)
- Alky Solidifier (15 oz. Bottle)
- Green-Z® Solidifier (15 oz. Bottle)
- Petro Lock Absorbent (8 oz. Solidifier)
- Caution Tape Roll
- 3 Clear Disposal Bags (30" x 37")\*
- D.O.T. Incident Report Tag
- Drop Cloth
- Instructions
- Pair of Nitrile Gloves
- Respirator
- 3 Scoops/Scrapers
- 3 Twist Ties
- Vented Goggles
- Zammie® Sewer Stopper (15" x 21")

# CONTAIN SPILLS WITH EASE

Meets requirements for compliance with DOT 1713 hazardous discharge; EPA 623.30 immediate action and 263.31 discharge cleanup

\*\*\* Careful consideration should be given to chemical spill cleanup, i.e., room size, temperature, quantity and concentration(s) of the solution(s).

\*Dispose of potentially infectious materials in accordance with all local, state and federal regulations.



Rev 25.08

 **Safetec**  
of America, Inc.

 **Made in the USA**

**Safetec of America, Inc.**  
1-800-456-7077 | [safetec.com](http://safetec.com)