



02/2024

## Corporate Statement

### RE: Shelf Stability and Expiration of SaniZide Plus® Germicidal Solution formula

Safetec's SaniZide Plus® Germicidal Solution formula will remain stable when stored in an approved container under the correct conditions and is not required to have expiration dating.

SaniZide Plus® Germicidal Solution formula is classified as a stable chemistry with both "accelerated and true time stability studies" that substantiate this formulation remains stable when stored in an approved container. Because it remains stable it does not display an expiration date.

#### ***In the "Storage and Disposal" section of the label on each container***

Do Not Contaminate Water, Food, Or Feed By Storage Or Disposal.

**PESTICIDE STORAGE** – Store In A Dry Place No Lower In Temperature Than 50°F Or Higher Than 120° F.

**PESTICIDE DISPOSAL** – Pesticide Wastes Are Acutely Hazardous. Improper Disposal Of Excess Pesticide, Spray Mixture, Or Rinsate Is A Violation Of Federal Law. If These Wastes Cannot Be Disposed Of By Use According To Label Instructions, Contact Your State Pesticide Or Environmental Control Agency, Or The Hazardous Waste Representative At The Nearest Epa Regional Office For Guidance.

**CONTAINER DISPOSAL**– Nonrefillable Container. Do Not Reuse Or Refill Container. Clean Container Promptly After Emptying. Triple Rinse As Follows: Fill Container ¼ Full With Water And Recap. Agitate Vigorously. Follow Pesticide Disposal Instructions For Rinsate Disposal. Drain For 10 Seconds After The Flow Begins To Drip. Repeat Procedure Two More Times. Offer For Recycling Or Reconditioning, If Available. If Not Available, Puncture And Dispose In A Sanitary Landfill.

#### **As per SECTION 7. HANDLING AND STORAGE in the SDS Sheet:**

##### ***Precautions for safe handling:***

*Handle in accordance with good industrial hygiene and safety practice. Do not ingest.  
Avoid contact with eyes.*

##### ***Conditions for safe storage :***

*Keep container tightly closed. Keep cool. Store away from incompatible materials.*

##### ***Incompatible materials :***

*Strong oxidizing agents. Acids.*

It is ultimately the responsibility of the customer to determine whether their inventory of SaniZide Plus® Germicidal Solution formula have been stored in accordance with the recommended conditions for safe storage and efficacy as outlined in our guidelines.



**02/2024**

## **Corporate Statement**

### **RE: Co-mixing (Topping off) SaniZide Plus® Germicidal Solution (EPA Reg. No. 1839-83-67161) Relative to Batch/Lot Numbers and Expiration Dates**

Clarification Regarding the Co-mixing (*Topping-off*) of SaniZide Plus® Germicidal Solution Relative to Batch Numbers and Expiration Dates.

The US EPA regulates the labeling and use of disinfectants used in the United States. Approved disinfectant product labels display an EPA Registration Number indicating that the product is approved for sale and use in the United States as well as EPA approval of the claims printed on the label. Also displayed on the container label is an Establishment Number that indicates the origin of the production facility as well as indicating the EPA's approval to produce the product at that facility. In addition, all of the fifty states require approval of disinfectants used in that state. It is a violation of federal law to use a product inconsistent with its labeled instructions.

In regards to the SaniZide Plus® Germicidal Solution: A batch/ lot number appears on the neck of each container or in some cases on the bottom of the container. The batch/ lot number indicates the specific product history including the date of production.

SaniZide Plus® Germicidal Solution is classified as a stable chemistry with both "accelerated and true time stability studies" that substantiate this product remains stable when stored in an approved container. Because it remains stable it does not display an expiration date.

An important note is that in the section labeled "Storage and Disposal" on each container, there are specific directions indicating as to whether the product container is approved for refilling or not. Unlike solutions that display expiration dates, and where the label specifically dictates that the container may not be refilled, the SaniZide Plus® Germicidal Solution container may be refilled from a larger container to a smaller dispenser container and the smaller container may be "top filled" thus co-mixing the solution.

When questioning whether a disinfectant may or may not be refilled into smaller containers always consider the expiration date and read the product Storage and Disposal directions.

*For additional clarification or questions contact Safetec of America, Monday through Friday at 1-800-456-7077.*