



+



Universal Spill Response



Multi-Use

Spill Leader Kit

Solidify various spills

This wall-mounted spill kit is designed for emergency response or hazardous lab spills. The solid solutions neutralize, stabilize, and solidify a variety of aqueous spills wherever they occur.

Acid Lock® Solidifier

Rapidly neutralizes and solidifies low pH spills:
Acetic Acid, Hydrochloric Acid, Nitric Acid, and Sulfuric Acid

Alky® Solidifier

Rapidly neutralizes and solidifies high pH spills:
Ammonium Hydroxide, Sodium Hydroxide, and Potassium Hydroxide

Green-Z® Solidifier

Solidifies blood and body fluids, as well as chemical*** spills

Petro Lock Solidifier

Rapidly absorbs petroleum-based spills:
Gasoline, Solvents, Crude Oil, and Paint Sludge

Red Z® Solidifier (contains chlorine)

Solidifies potentially infectious blood and body fluids

What's Included:

- Adhesive Wall Mount
- Acid Lock® Solidifier (5 oz. Bottle)
- Alky® Solidifier (5 oz. Bottle)
- Green-Z® Solidifier (5 oz. Bottle)
- Petro Lock Absorbent (2 oz. Solidifier)
- Red Z® Solidifier (5 oz. Bottle)
- D.O.T. Incident Report Tag
- 2 Hazardous Waste D.O.T. Labels
- 3 Pairs of Nitrile Gloves
- 2 Non-Hazardous Waste D.O.T. Labels
- 6 Scoops/Scrapers
- 6 Zip Lock Bags (9" x 12")

Order No.	Container	Container Dimensions	Case Qty.	Case Wt.
15205	Poly Bag	14.125" x 3.125" x 20.5"	1/case	5 lbs.

THE ULTIMATE SPILL CONTROL SOLUTION

**Meets requirements for compliance with DOT
1713 hazardous discharge; EPA 623.30 immediate
action and 263.31 discharge cleanup**

*** Careful consideration should be given to chemical spill cleanup, i.e., room size, temperature, quantity and concentration(s) of the solution(s).

*Dispose of potentially infectious materials in accordance with all local, state and federal regulations.



Safetec
of America, Inc.

Made in the USA

Safetec of America, Inc.
1-800-456-7077 | safetec.com