



+



Universal Spill Response



Multi-Use

## Spill Leader Kit

### Solidify Various Spills

This wall-mounted spill kit is designed for emergency response and hazardous lab spills. The solidifier solutions quickly neutralize, stabilize, and solidify a variety of aqueous spills such as chemicals, body fluids, or other hazardous substances, wherever they occur.

#### Acid Lock Solidifier

Rapidly neutralizes and solidifies **low pH spills**:  
Acetic Acid, Hydrochloric Acid, Nitric Acid, and Sulfuric Acid

#### Alky Solidifier

Rapidly neutralizes and solidifies **high pH spills**:  
Ammonium Hydroxide, Sodium Hydroxide, and Potassium Hydroxide

#### Green-Z® Solidifier

Solidifies **blood and body fluids**, as well as **chemical\*\*\* spills**

#### Petro Lock Solidifier

Rapidly absorbs **petroleum-based spills**:  
Gasoline, Solvents, Crude Oil, and Paint Sludge

#### Red Z® Solidifier (contains chlorine)

Solidifies potentially infectious **blood and body fluids**

## What's Included:

- Adhesive Wall Mount
- Acid Lock Solidifier (5 oz. Bottle)
- Alky Solidifier (5 oz. Bottle)
- Green-Z® Solidifier (5 oz. Bottle)
- Petro Lock Absorbent (2 oz. Solidifier)
- Red Z® Solidifier (5 oz. Bottle)
- D.O.T. Incident Report Tag
- 2 Hazardous Waste D.O.T. Labels
- 3 Pairs of Nitrile Gloves
- 2 Non-Hazardous Waste D.O.T. Labels
- 6 Scoops/Scrapers
- 6 Zip Lock Bags (9" x 12")

Order No.	Container	Container Dimensions	Case Qty.	Case Wt.
15205	Plastic Wall Mount	14.125" x 3.125" x 20.5"	1/case	5 lbs.

# THE ULTIMATE SPILL CONTROL SOLUTION

Meets requirements for compliance with DOT 1713 hazardous discharge; EPA 623.30 immediate action and 263.31 discharge cleanup

\*\*\* Careful consideration should be given to chemical spill cleanup, i.e., room size, temperature, quantity and concentration(s) of the solution(s).

\*Dispose of potentially infectious materials in accordance with all local, state and federal regulations.



Rev 25.08

 **Safetec**  
of America, Inc.

 **Made in the USA**

**Safetec of America, Inc.**  
1-800-456-7077 | safetec.com