



+



Universal Spill Response



Multi-Use

Spill Leader Kit

Solidify various spills

This wall-mounted spill kit is designed for emergency response or hazardous lab spills. The solid solutions neutralize, stabilize, and solidify a variety of aqueous spills wherever they occur.

Acid Lock® Solidifier

Rapidly neutralizes and solidifies low pH spills:

Acetic acid, Hydrochloric acid, Nitric acid, and Sulfuric acid

Alky® Solidifier

Rapidly neutralizes and solidifies high pH spills:

Ammonium Hydroxide, Sodium Hydroxide, and Potassium Hydroxide

Green-Z® Solidifier

Solidifies blood and body fluids, as well as chemical*** spills

Petro Lock Solidifier

Rapidly absorbs petroleum-based spills:

Gasoline, Solvents, Crude Oil, and Paint Sludge

Red Z® Solidifier (contains chlorine)

Solidifies potentially infectious blood, and body fluids

Order No.	Description	Case Qty.	Case Wt.
15205	Spill Leader Kit	1/case	5 lbs.

What's included:

- Acid Lock® Solidifier (5 oz. Bottle)
- Alky® Solidifier (5 oz. Bottle)
- Green-Z® Solidifier (5 oz. Bottle)
- Petro Lock Absorbent (2 oz. Solidifier)
- Red Z® Solidifier (5 oz. Bottle)
- D.O.T. Incident Report Tag
- Hazardous Waste D.O.T. Label (2)
- Nitrile Gloves (3)
- Non-Hazardous Waste D.O.T. Label (2)
- Respirator
- Scoop/Scraper (6)
- Zip Lock Bags (6 - 9" x 12")

Unopened Hard Case Dimensions: 3.375" x 14.375" x 20.5"

THE ULTIMATE SPILL CONTROL SOLUTION

Meets requirements for compliance with DOT 1713 hazardous discharge; EPA 623.30 immediate action and 263.31 discharge cleanup.

*** Careful consideration should be given to chemical spill cleanup, i.e., room size, temperature, quantity and concentration(s) of the solution(s).

*Dispose of potentially infectious materials in accordance with all local, state and federal regulations.

Safetec
of America, Inc.



Made in the USA

Safetec of America, Inc.
1-800-456-7077 | safetec.com