



Drainage

Available in:



Multiple Sizes

Continent Catheters

Torbot® Continent Catheters are designed for dependable drainage of internal continent reservoirs created through continent ostomy and continent urostomy procedures. Available in 24 French and 30 French sizes, these flexible straight-tip catheters support routine intermittent catheterization while accommodating different drainage needs and surgical procedures. The soft, straight-tip design provides a comfortable alternative to rigid bullet-tip catheters, helping support everyday patient care.

Directions For Use:

Cut the catheter into two pieces using Ostomy Scissors. Connect the sections using the coupling piece and empty the reservoir. Once emptied, disconnect the connector and clean the catheter sections. Store back into the plastic bag.

Ideal for use in:

- BCIR Catheterization
- Continent Ostomy Drainage
- Continent Urostomy Drainage
- Indiana Pouch Catheterization
- Kock (K-Pouch) Catheterization
- Mitrofanoff Catheterization
- Routine Intermittent Catheterization Through a Stoma

Order No.	Size	Product Type	Bag	Case Qty.	Case Wt.
ME8724	24 Fr	Continent	25 ct.	6/case	8 lbs.
ME8724-2PK	24 Fr	Continent	25 ct.	2/case	3 lbs.
ME8730	30 Fr	Continent	25 ct.	6/case	8 lbs.
ME8730-2PK	30 Fr	Continent	25 ct.	2/case	3 lbs.

Catheter selection should always follow the healthcare provider's recommendation.

DESIGNED FOR CONTINENT DIVERSIONS

Flexible Straight-Tip Design

Available in Two Sizes



Rev 26.08