



Urinary Night Drainage Kit

The Torbot® Urinary Night Drainage Kit is designed to provide dependable overnight urinary drainage for patients with continuous urinary drainage needs. The complete kit includes a wide-base collection bottle, 5-foot drainage hose with connector, and one bottle of T3 Odor Eliminator. Designed with a lower center of gravity to help reduce tipping, it provides stable bedside urine collection for standard urostomy patients with an external urinary pouch and patients using Foley or male external catheter drainage systems.

Ideal for use in:

- Bedside Urine Collection
- Home Healthcare
- Long-Term Urinary Drainage Management
- Overnight Urinary Drainage
- Standard Urostomy Drainage Support
- Urinary Catheter Drainage Support

Order No.	Container	Container Dimensions	Case Qty.	Case Wt.
TSN4326	Plastic Bag	1" x 12" x 15"	4/case	3 lbs.

STABLE OVERNIGHT DRAINAGE

Wide-Base Stable Bottle

Designed for Overnight Drainage



Rev 26.08