



# Universal Precaution Compliance Kit

This kit is crucial due to the health risks associated with bloodborne pathogens like HBV, HCV, and HIV. It adheres to **OSHA, CDC, and State Health Department** standards, providing the necessary personal protective equipment (PPE) and cleanup tools.

Under the OSHA compliance standard, when there is occupational exposure to blood and body fluids, employers are required to provide appropriate, readily accessible PPE at no cost to employees. Contaminated work surfaces must be decontaminated with an appropriate disinfectant after a spill cleanup. Additionally, a means to sanitize hands and other skin must be available.

Suitable for various workplaces, such as nursing homes, hospitals, schools, correctional facilities, and any business where job-related contact with blood is a potential risk.

This kit includes Red Z<sup>®</sup> Solidifier, a quick and efficient solution for solidifying blood and body fluids while eliminating odors.

## What's Included:

- 10 g Red Z<sup>®</sup> Solidifier
- Combo Mask/Safety Shield
- Identification Tag
- Instructions
- Pair of Vinyl Gloves
- p.a.w.s.<sup>®</sup> Antimicrobial Hand Wipe
- Protective Apron
- Red Biohazard Bag (24" x 24")\*
- SaniZide Pro 1<sup>®</sup> Surface Disinfectant Wipe
- Scoop/Scraper
- Twist Tie

Order No.	Container	Container Dimensions	Case Qty.	Case Wt.
17100	Poly Bag	1" x 6.19" x 9.75"	24/case	6 lbs.
17102	Hard Case	2.63" x 6.5" x 8.12"	12/case	10 lbs.

# ACCIDENTS HAPPEN

Meets OSHA Compliance 1910.1030 – Bloodborne Pathogens Standard



\*Dispose of potentially infectious materials in accordance with all local, state and federal regulations.

Rev 25.08



Safetec of America, Inc.  
1-800-456-7077 | safetec.com