



**Fast-Acting**



**Ready-To-Use**



**Multiple Sizes**

## Zorb Sheets

### *Blood, Body, & Chemical Fluid Absorbents*

Transfer specimens safely and easily using Safetec's Zorb Sheets! These versatile sheets help reduce compliance risks associated with the handling, shipping, diagnostic specimens, infectious substances, hazardous materials and dangerous goods. The tissue-based material allows fluid to immediately activate the superabsorbent polymers to solidify and transform the contents into a semi-solid (gel) mass. Safetec's Zorb Sheets—your reliable choice for safe transportation that helps comply with DOT, CDC, OSHA, IATA, ICAO, USPS, and UN regulations.†

- Available in multiple sizes
- Fast-acting
- Easy-to-use
- Lightweight
- Space saving
- Prevents leaking

#### **Ideal for use in:**

- **Healthcare**
  - Biohazard Bags
  - Disposable Lab Containers
  - Glassware
  - Patient Testing Centers
  - Sharps Containers
  - Triple Packaging Systems
- **Medical Device Manufacturers**
- **Medical Shippers**
- **Occupational Health & Safety**
  - Crime Scene Investigation
  - Mortuary Industry
  - Shipping Facilities  
(Packaging regulations)
  - Tattoo Shops  
(Ink containers)

# COMBAT SPILLS

**For diagnostic specimens, infectious substances, and dangerous goods**

**Zorb Sheets help comply with DOT, CDC, OSHA, IATA, ICAO, USPS and UN regulations†**

♦ Data is an approximation of values. Blood and urine solidification are based on the industry standard 0.9% NaCl saline solution.

\* Dispose of infectious and hazardous waste in accordance with all local, state and federal regulations.

† Visit our website for a full list of regulation references.



**Made in the USA**

**Safetec of America, Inc.**  
1-800-456-7077 | [safetec.com](https://www.safetec.com)



# Zorb Sheet Types

## S/S Zorb Sheets

*Single-Ply with Standard Absorbency*

The super absorbent polymer in Single Laid (S/S) sheets is pressed into a single sheet, providing efficient absorbency.




Order No.	Size	Max. Blood/Urine/Saline Absorbency*	Max. Water Absorbency*	Bags	Case Qty.	Case Wt.
44001	3" x 3"	(Up to 20 cc)	(Up to 131 cc)	1,000 ct.	20/case	28 lbs.
44002	12" x 12"	(Up to 320 cc)	(Up to 2,096 cc)	-	500/case	11 lbs.
44003	6" x 6"	(Up to 80 cc)	(Up to 524 cc)	500 ct.	5/case	13 lbs.
44004	3" x 6"	(Up to 40 cc)	(Up to 262 cc)	500 ct.	5/case	7 lbs.

## A/S Zorb Sheets

*Multi-Ply with Extra Absorbency*

Air Laid (A/S) sheets feature the super absorbent polymer pressed between two sheets, offering enhanced absorbency and even more effective containment.



Order No.	Size	Max. Blood/Urine/Saline Absorbency*	Max. Water Absorbency*	Bags	Case Qty.	Case Wt.
44016	12" x 12"	(Up to 480 cc)	(Up to 2,240 cc)	-	500/case	18 lbs.
44020	2.5" x 4.25"	(Up to 35 cc)	(Up to 165 cc)	100 ct.	50/case	14 lbs.
44035	5" x 6"	(Up to 100 cc)	(Up to 467 cc)	500 ct.	10/case	37 lbs.
44062	4" x 5"	(Up to 67 cc)	(Up to 311 cc)	100 ct.	10/case	5 lbs.
44090	6" x 12"	(Up to 240 cc)	(Up to 1,120 cc)	-	1,000/case	10 lbs.
45002	4" x 3"	(Up to 40 cc)	(Up to 187 cc)	1,000 ct.	10/case	31 lbs.
45009	9" x 12"	(Up to 360 cc)	(Up to 1,680 cc)	-	500/case	13 lbs.
45010	8" x 9"	(Up to 240 cc)	(Up to 1,120 cc)	125 ct.	8/case	17 lbs.

Custom sizes and imprints are available. Call for a quote.

\*Data is an approximation of values. Blood and urine solidification are based on the industry standard 0.9% NaCl saline solution.

†Visit our website for a full list of regulation references.