



Fast-Acting



Ready-To-Use



Multiple Sizes

Zorb Sheets

(Blood/Body & Chemical Fluid Absorbents)

The tissue-based material allows fluid to immediately activate the superabsorbent polymers to solidify and transform the contents into a semi-solid (gel) mass. In the event of a spill or leak, it creates an extra protection barrier for workers' safety. (These high quality granules have the unique ability to hold large amounts of absorbed water-based fluid, even against pressure.)

- Available in multiple sizes (custom sizes available)
- Fast-acting and easy-to-use
- Lightweight (space saving)
- Prevents leaking

Ideal for use in:

- **Healthcare**
 - Biohazard Bags
 - Disposable Lab Containers
 - Glassware
 - Patient Testing Centers
 - Sharps Containers
 - Triple Packaging Systems
- **Medical Device Manufacturers**
- **Medical Shippers**
- **Occupational Health & Safety**
 - Crime Scene Investigation
 - Mortuary Industry
 - Shipping Facilities (Packaging regulations)
 - Tattoo Shops (Ink Containers)

Order No.	Description	Sheet Size	Case Qty.	Case Wt.	Max. Blood/Urine/Saline Absorbency [†]	Max. Water Absorbency [†]
44001	Zorb Sheet	3" x 3"	20,000	26 lbs.	(Up to 20cc)	(Up to 60cc)
44002	Zorb Sheet	12" x 12"	500	11 lbs.	(Up to 275cc)	(Up to 720cc)
44003	Zorb Sheet	6" x 6"	2,500	14 lbs.	(Up to 80cc)	(Up to 240cc)
44004	Zorb Sheet	3" x 6"	2,500	3.5 lbs.	(Up to 36cc)	(Up to 120cc)

COMBAT SPILLS

For diagnostic specimens, infectious substances, and dangerous goods

Zorb Sheets help comply with DOT, CDC, OSHA, IATA, ICAO, USPS and UN regulations†

[†]Data is an approximation of values. Blood and urine solidification are based on the industry standard 0.9% NaCl saline solution.

* Dispose of infectious and hazardous waste in accordance with all local, state and federal regulations.

†Visit our website for a full list of regulation references.



Safetec of America, Inc.
1-800-456-7077 | [safetec.com](https://www.safetec.com)