



"Cleans, disinfects and deodorizes"



- Kills Pandemic 2009 H1N1 Influenza A
- Effective Against TB, MRSA, HIV-1, Norovirus and More
- EPA Registered
- Exceptionally Durable, Pure and Low Linting Wipes
- Alcohol Free
- Cleans, Disinfects and Deodorizes
- For Use on All Hard, Non-Porous Surfaces Normally Treated with a Disinfectant*

**For plastic and painted surfaces, spot test on an inconspicuous area before use.*

Safetec® Sanizide Plus® Germicidal Wipes

Sanizide Plus® Germicidal Wipes are pre-saturated hard surface wipes that kill TB, MRSA, VRE, E. Coli and more in 10 minutes and HIV-1 in 2 minutes. When you reach for a Sanizide Plus® Wipe you get a wipe that is clean wipe after wipe. The towelettes leave behind virtually no lint or particles to scratch or contaminate surfaces. Our dual quat, alcohol free formula is non-flammable, safe to ship and non-corrosive on surfaces. Sanizide Plus® Wipes are ready-to-use, convenient and user-friendly.

Order No.	Description	Case Qty.
34823	160ct. Sanizide Plus® Germicidal Wipe Canister	12 per cs.
34826	50ct. Sanizide Plus® Germicidal Wipe Dispenser Box	6 boxes per cs.
2510097	160ct. Canister Holder	1 each

