



Chemotherapy Spill Kit

Protect and Clean up Spills Involving Cytotoxic Therapy Drugs. This kit features Green-Z[®] solidifier to solidify and contain hazardous chemical spills for quick and easy clean-up and disposal with minimal exposure.

Item#	Description	Case Qty.	Case Wt.
49725	Chemotherapy Spill Kit (Poly Bag)	12	8 lbs.

Compliance

*Nitrile gloves tested for use with Chemotherapy drugs using ASTM D6978. (They are powder free, latex-free, DEHP free and silicone free)

**The Chemotherapy gown is compatible with OSHA recommendations, and ASHP and ONS guidelines for wearing a gown while preparing or administering chemotherapy drugs. The gown is disposable, lint-free, made of low permeability fabric, constructed with a solid front and back closure and is long sleeved with tight-fitting cuffs.



Each Chemotherapy Spill Kit Contains:

- 2 Pairs of Nitrile Chemotherapy Gloves*
- Full-Face Safety Shield
- Chemotherapy Gown**
- Green-Z[®] Solidifier (21g pouch)
- Scoop & Scraper
- Wiper Pads (2)
- Shoe Covers (non-slip)
- Red Bio-hazard Waste Bag (33g)
- Twist Tie
- Identification Tag
- Instruction Card
- Caution Sign