



Formaldehyde Spill Response Kit & FSC-1

Reduce Exposure to Formaldehyde Gas

This kit is a convenient emergency response kit for medical personnel using formalin or other formaldehyde containing substances. It is ideal for containing small volume spills. (Where a significant volume of fluid must be contained, notify emergency response personnel in accordance with facility guidelines.) Available in highly visible, durable, easy-to-share hard vinyl case.

Item#	Description	Case Qty.	Case Wt.
48630	Formaldehyde Spill Kit (Hard Case)	1	5 lbs.
48625	Formaldehyde Spill Kit (Kit Refill)	1	4 lbs.
47101	11 oz. FSC-1 Shaker Bottle	12	11 lbs.



Each Formaldehyde Spill Response Kit Contains:

- FSC-1 Solidifier (11 oz bottle)
- Nitrile Gloves
- Vented Goggles
- Scoop & Scraper (4)
- Respirator
- Wiper Pads (3)
- Open Backed Gown
- Gallon Bags (4 – 30" x 36")



FSC-1 Formalin Spill Control

This aids in reducing the exposure to formaldehyde gas by immobilizing formalin and transforming it to an easily removed semi-solid mass. The resulting mass cools rapidly to minimize formaldehyde off-gassing.

*FSC-1 capacity to convert and immobilize formalin decreases as the concentration of formalin increases.

Formaldehyde Strength (% wt.)	Absorption Capacity (per 11oz. jar)	
	milliliters (mL)	fluid ounces (fl. oz.)
4%	1800	60.8
10%	720	24.3
15%	480	16.2
20%	360	12.1
30%	240	8.1
37%	195	6.5

Please Note: Not recommended to be used with 40% formalin or higher.