



Is Your First Aid Kit ANSI Compliant?

Class A First Aid Kit Checklist

First Aid Kit Supply	Minimum Quantity	Minimum Size/Vol.	Do You Have It?
Adhesive Bandage	16	1 x 3 in.	
Adhesive Tape	1	2.5 yd (total)	
Antibiotic Application 	10	1/57 oz. (0.5 g)	
Antiseptic 	10	1/57 oz. (0.5 g)	
Breathing Barrier	1	N/A	
Burn Dressing (gel soaked)	1	4 x 4 in.	
Burn Treatment 	10	1/32 oz. (0.9 g)	
Cold Pack	1	4 x 5 in.	
Eye Covering (with means of attachment)	2	2.9 sq. in.	
Eye/Skin Wash	1	1 fl oz. total	
Foil Blanket	1	52 x 84 in.	
First Aid Guide	1	N/A	
Hand Sanitizer 	10	1/32 oz. (0.9 g)	
Medical Exam Gloves	2 pairs	N/A	
Roller Bandage	1	2 in. x 4 yd	
Scissors	1	N/A	
Splint	0	4 x 24 in.	
Sterile Pad	2	3 x 3 in.	
Tourniquet	0	1.5 in. (width)	
Trauma Pad	2	5 x 9 in.	
Triangular Bandage	1	40 x 40 x 56 in.	

*Information based on the [ANSI/ISEA Z308.1-2021](#) American National Standards for Minimum Requirements for Workplace First Aid Kits and Supplies.

 Safetec product ready to ship



Missing something? We've got you covered.

ANSI compliant Safetec products ready to ship:

- Triple Antibiotic and Single Antibiotic Ointments
- Antiseptic Wipes and Spray
- First Aid & Burn Cream, Burn Gel, and Burn Spray
- Instant Hand Sanitizer







Note: Safetec does not manufacture pre-packaged first aid kits. However, we can assemble kits using Safetec supplies and/or customer-provided supplies.



Safetec of America, Inc.
1-800-456-7077 | safetec.com

Is Your First Aid Kit ANSI Compliant?

Class B First Aid Kit Checklist

First Aid Kit Supply	Minimum Quantity	Minimum Size/Vol.	Do You Have It?
Adhesive Bandage	50	1 x 3 in.	
Adhesive Tape	2	2.5 yd (total)	
Antibiotic Application 	25	1/57 oz. (0.5 g)	
Antiseptic 	50	1/57 oz. (0.5 g)	
Breathing Barrier	1	N/A	
Burn Dressing (gel soaked)	2	4 x 4 in.	
Burn Treatment 	25	1/32 oz. (0.9 g)	
Cold Pack	2	4 x 5 in.	
Eye Covering (with means of attachment)	2	2.9 sq. in.	
Eye/Skin Wash	1	4 fl oz. total	
Foil Blanket	1	52 x 84 in.	
First Aid Guide	1	N/A	
Hand Sanitizer 	20	1/32 oz. (0.9 g)	
Medical Exam Gloves	4 pairs	N/A	
Roller Bandage	2	2 in. x 4 yd	
	1	4 in. x 4 yd	
Scissors	1	N/A	
Splint	1	4 x 24 in.	
Sterile Pad	4	3 x 3 in.	
Tourniquet	1	1.5 in. (width)	
Trauma Pad	4	5 x 9 in.	
Triangular Bandage	2	40 x 40 x 56 in.	

*Information based on the [ANSI/ISEA Z308.1-2021](#) American National Standards for Minimum Requirements for Workplace First Aid Kits and Supplies.

 Safetec product ready to ship



Missing something? We've got you covered.

ANSI compliant Safetec products ready to ship:

- Triple Antibiotic and Single Antibiotic Ointments
- Antiseptic Wipes and Spray
- First Aid & Burn Cream, Burn Gel, and Burn Spray
- Instant Hand Sanitizer



Note: Safetec does not manufacture pre-packaged first aid kits. However, we can assemble kits using Safetec supplies and/or customer-provided supplies.

Safetec of America, Inc.
1-800-456-7077 | safetec.com